

interscience



Anabac[®]

Autoclave deodorant

Anabac[®] Brochure N° 30910-0



Give some air to your lab!

Laboratories are very concerned about hygiene cleanliness and comfort of work. Try **Anabac®**, autoclave deodorant, for a perfumed touch in your lab. **Anabac®** autoclave deodorant bring the discrete final touch, the refinement of cleanliness.



Classic

Well perfumed
All applications

• *mint and eucalyptus fragrance*



Poma

Sweet perfume
Food industry

• *apple fragrance*



Citrus

Sweet perfume
Food industry

• *lemon fragrance*



Floral

Well perfumed
All applications

• *flower fragrance*




Peach

Well perfumed
All applications

• *peach fragrance*




1.




During sample sterilization, bad smell can emerge from the autoclave.

2.



Place one capsule (or 2 if the odor is strong) in the water at the bottom of the autoclave.

3.



After a sterilization cycle, the capsule melts, allowing the perfume to be released.

Advantages

- *Work comfort*
- *Compatible with any autoclave*
- *No effect on the autoclaving*
- *Pleasant perfume*
- *100% biodegradable*

Applications

- For any kind of autoclave (max: 200 L), recommended for food industry.

Specifications

- Pot of 100 capsules
- Pot weight: 160 g
- 5 fragrances: - Classic (ref: 320 100)
- Poma (ref: 320 200)
- Citrus (ref: 320 300)
- Floral (ref: 320 400)
- Peach (ref: 320 500)



Please return this form to your local distributor or to interscience.



In addition

Anabac® autoclave deodorant

- Quotation
- Free "5 fragrances" sample

I would like to receive more info on:

- BagMixer® 400 ml** blender
- JumboMix® 3500 ml** blender
- BagFilter®** filter bags
- BagTools®** accessories
- Gravimat®** gravimetric dilutor
- easySpiral®** automatic plater
- Scan®** colony counters

Name.....

Company.....

Address.....

.....

Zip/City..... Country.....

Phone..... Fax.....

email.....

web.....



www.ert.pt
info@ert.pt

Telf. 217 142 066 - Fax. 217 145 674

Rua Soeiro Pereira Gomes, Lote 5, Cave Drt.
2845-387 Amora Portugal

interscience

Worldwide: 30. Ch Bois Arpents, 78860 St Nom FRANCE
Tel. +33(0)1 34 62 62 61 / Fax. +33(0)1 34 62 43 03
info@interscience.fr / www.interscience.fr

USA/Canada: 199 Weymouth St. ROCKLAND, MA 02370 USA
Tel. +1 781 792 2133 / Fax. +1 781 792 2134
info@intersciencelab.com / www.intersciencelab.com



Clinet® 101

Chlorhexidine digluconate disinfecting wipes

Ideal for laboratories, hospitals, samplers, home use.



Clinet® 202

Individually packed chlorhexidine digluconate disinfecting wipes

Ideal for laboratories, hospitals, samplers, transports, hand hygiene.



Clinet® 303

Bactericidal and fungicidal spray

Ideal for large surfaces, hard to reach places, angle of walls, rough edges.



Clinet® 404 SaniFood®

Food contact approved disinfecting wipes

For food professionals: disinfection without rinsing of surfaces and materials in contact with food.



Clinet® Anti-algae

Algicide for laboratories water baths.

1 liter bottle with incorporated doser, Clinet®Anti-algae is used diluted at 0.1% to inhibit the growth of algae in water baths.

Anabac® Brochure N° 30910-0 - Anabac®, Clinet®, Sanifood®, BagClip®, BagMixer®, JumboMix®, BagFilter®, BagTools®, easySpiral® and Scan® are registered trademarks of interscience. Information in this leaflet is not contractually binding.