








So Easy, So Charming




-  **Focus on nucleic acid:** 230nm/260nm/280nm
-  **Convenient:** Only 2 μ l sample needed
-  **Accurate:** High reproducibility and repeatability
-  **All-in-One:** Embedded computer and thermo printer
-  **Economic:** Much less cost for most useful functions



MAESTRO nano

Specifications

Dimension	275X275X110 mm
Concentration	2~2000 ng/μl
Sample Size	2 μl
Voltage	100~240V
Light	Xenon Lamp
Data Output	Printer/USB



www.ert.pt

info@ert.pt

Telf. 217 142 066 - Fax. 217 145 674

Rua Soeiro Pereira Gomes, Lote 5, Cave Drt.
2845-387 Amora Portugal

MAESTROGEN

www.maestrogen.com

info@maestrogen.com

