

Clinet® 202

Individually-packed cleaning wipes

THE FINISHING TOUCH

Individually-packed Clinet® 202 cleaning wipes can be used wherever a high level of hygiene is needed. Their compact format allows them to be carried in a pocket, bag, car glove compartment... and are immediately available in any emergency.


CHARACTERISTICS

The disinfectant used is an alcoholic solution of benzalchonium chloride and chlorhexidine digluconate.

Renowned for its remarkable effectiveness, Clinet® 202 solution is hypoallergenic and complies to the following norms :

- Bactericidal : NF EN 1040 - NFT T 72.151 - NF T 72.170 - NFT T 72.190
- Fungicidal : NF EN 1275 - NF T 72.201 - NF T 72.190
- Virucidal : Active on Herpes, Hepatitis B and HIV 1 (in vitro).

Clinet® 202 are ready for use. Simply wipe down the surface to be disinfected. No need to rinse.



Individually-packed wipes.
Wipe size : 20 x 30 cm.
Box of 100 wipes.
1.4 kg - 21.5 x 26 x 10.5 cm.

Precautions :
For external use only.
Do not apply on mucous membranes or baby skin.

For cleaning and disinfecting large surfaces in the laboratory, production shop and all places of public affluence.

Production line, laboratory equipment, lavatory, telephone, ambulance, car, stretcher, public places, hands...

.....

Please send me : a price offer for Clinet® 202.
 technical information on Clinet® 202.
 a free sample of Clinet® 202.

Company

Name

Position

Address

Zip code/City

Country

Tel. Fax

E-mail

Hygiene & Security

- | | info. | samples |
|---|--------------------------|--------------------------|
| Clinet® 101 (chlorhexidine cleaning wipes) | <input type="checkbox"/> | <input type="checkbox"/> |
| Clinet® 303 (bactericidal spray) | <input type="checkbox"/> | <input type="checkbox"/> |
| Clinet® anti-algae (algae inhibiting liquid) | <input type="checkbox"/> | <input type="checkbox"/> |
| Anabac® (autoclave deodorant) | <input type="checkbox"/> | <input type="checkbox"/> |
| BioParfum® (odor neutralizers) | <input type="checkbox"/> | <input type="checkbox"/> |

BagSystem® for sample preparation

- | | | |
|--|--------------------------|--------------------------|
| BagMixer® (sample blenders) | <input type="checkbox"/> | |
| BagFilter® (sterile blender bags) | <input type="checkbox"/> | <input type="checkbox"/> |
| BagPipet® (automatic pipette) | <input type="checkbox"/> | |

www.ert.pt
info@ert.pt

Tel. 217 142 066 - Fax. 217 145 674

Rua Soeiro Pereira Gomes, Lote 5, Cave Drt.
2845-387 Amora Portugal

Please fax this request to your local dealer or to :